# eldom sonic

SD200

SONICZNA SZCZOTECZKA /**PL**/  
SONIC TOOTHBRUSH /**EN**/  
SONICKÁ ZUBNÍ KARTÁČEK /**CZ**/  
SONIC ZAHRNBÜRSTE /**DE**/  
ЗВУКОВАЯ ЗУБНАЯ ЩЕТКА /**RU**/  
SONICKÁ ZUBNÁ KEFKA /**SK**/  
SZŐNIKUS FOGKEFE /**HU**/  
SÓNICA CEPILLO DE DIENTES /**ES**/

**Eldom Sp. z o.o.** • ul. Pawła Chromika 5a • 40-238 Katowice, POLAND  
tel: +48 32 2553340 • fax: +48 32 2530412 • [www.eldom.eu](http://www.eldom.eu)

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**/EN/** Disposal of used electrical and electronic equipment (applies to the European Union countries and other European countries with separate waste collection systems).

The symbol on the product or its packaging indicates that the product should not be classified as household waste. It should be handed over to a collection point that deals with the collection and recycling of electrical and electronic devices. Correct disposal of the product prevents potential negative consequences for the environment and human health that could be caused by improper disposal of the product. Material recycling helps to preserve natural resources. For detailed information on recycling of this product, please contact your local authority, a cleaning service provider or a shop where the product has been bought.

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## **SAFETY INSTRUCTIONS**

- Before using the appliance for the first time, please read carefully these Instructions for Use.
- The device can be used by children aged 8 and older and by persons with limited physical, sensory or mental capabilities, or persons without sufficient experience and knowledge only under supervision or if previously instructed on the safe use of the appliance and the possible risks. Children shall not play with the appliance. Cleaning and maintenance that are the responsibility of the user shall not be carried out by children without supervision.
- Do not place the appliance in the vicinity of heat sources.
- The appliance is intended for domestic use only.
- Do not immerse the cord or the appliance in water.
- The appliance must always be cleaned after use. The appliance must always be turned off before cleaning and maintenance as well as in case of malfunction and after operation.
- Do not use with a broken or otherwise damaged hairdressing head as it may lead to injury.
- Do not use in the open air.
- Use only with the originally included accessories and feeder.
- Never use the appliance when it is wet, with wet hands or in containers.
- The appliance must not be immersed in or sprinkled with water.
- When in operation, do not leave the appliance unattended.
- The appliance can be repaired only in an authorised Service Centre. Any modernisation or application of spare parts or elements other than the original ones is forbidden and dangerous for the user.
- The company Eldom Sp. z o.o. shall not be liable for any damage caused by inappropriate use of the appliance.

### **NOTE:**

- This device is designed exclusively for cleaning teeth, gums and tongue.
- When using this product for the first time, some people may experience discomfort. Use the toothbrush normally and any unpleasant sensations should gradually disappear within a week without hurting your teeth or gums.
- If you have received dental treatment in the last 2 months, consult your dentist before using the toothbrush.
- If you experience bleeding after using the toothbrush, consult your dentist.

**GENERAL INFORMATION**

1. Spindle
2. Housing
3. Power switch
4. Operating mode indication
5. Empty battery indication
6. Cover
7. Replaceable head
8. Charging station (USB)

**INTENDED USE**

The SD200 toothbrush has been designed to use sonic technology. During operation, the toothbrush head vibrates at speeds ranging from 31,000 to as many as 42,000 beats per minute, dynamically cleaning teeth and inter-dental spaces. This makes it possible to remove impurities from hard-to-reach areas more accurately.

**CHARGING THE DEVICE**

**The device may only be used after it has been charged.**

**Switch off the device before placing it on the charging station (8).**

- Connect the charging station (8) to the power supply unit (not included) with the parameters provided above, and then to the mains.
- Place the device on the charging station (8).
- When charging, the light (5) will be on and will go off when the device is charged.
- The device cannot be used while charging.
- The toothbrush will switch off automatically after placing it on the charging station.

**TECHNICAL DATA**

**Power:** 3W

**Speed:** 31,000 - 42,000 bpm

**Battery:** DC 3.7V, 700mAh, Li-Ion

**Charging time:** 16 hours

**Operating time:** up to 20 days (2 minutes twice/day)

**Toothbrush water-tightness:** IPX7

**Material**

- housing: ABS/TPE

- head bristles: Dupont Nylon

**Charging method:** induction

Connect the station (8) to a power supply unit with the following parameters: DC 5V, 0.6mA (not included)

**REPLACING THE HEAD (7)**

The brush head (7) is a working part and should be replaced every 3 months for hygienic reasons and whenever its bristles are worn.

It is not subject to warranty repair. The brush head (7) can only be replaced by the original brush head type included in the SDK200 range.

- The head can only be replaced when the toothbrush is switched off.
- Replacing the head (7) while the toothbrush is running may cause damage to the spindle (1).
- Attach the head (7) as shown in Fig. A.
- Keep a small gap between the housing (2) and the head (7) to allow proper device operation. If there is no gap, the device will be damaged when turned on

**OPERATING MODE SELECTION**

- To select the operating mode, press the button (3) until the correct light (4) comes on.
- Press and hold the button (3) for 3 seconds to change the mode while the device is running.
- Press and hold the button (3) for more than 3 seconds to turn the device off.

- The toothbrush will always begin working in the last mode selected.

## USE

- Place the head (7) on the spindle (1) – Fig. A
- Switch on the device, select the operating mode (see above) and switch it off.
- Wet the bristles of the head (7).
- Apply a small amount of toothpaste to the wet bristles.
- Place the toothbrush in your mouth, set it at the right angle (Fig.C) and turn it on.
- Brush your teeth in circular motions as described below without applying too much pressure to the teeth and gums.
- The brushing cycle takes 2 minutes. The toothbrush will generate an impulse every 30 seconds.
- The device will switch off automatically at the end of the set time period.
- The device operation can be stopped at any moment by pressing the switch (3).
- Rinse the mouth after brushing. Specialist rinses may be used.
- Remove the head (7), wash it thoroughly and set aside to dry.
- Wipe the device and set it down in a dry place or on the charging station (8).

## TEETH CLEANING – RECOMMENDATIONS

- Divide the teeth into 4 areas – Fig. B.
- Brush each part for approximately 30 seconds (approximately 45 seconds in deep cleaning mode). The toothbrush signals when the time is up.
- Apply the head bristles to the teeth correctly – Fig. C

## CLEANING AND MAINTENANCE

Clean the appliance after each use.

- Before cleaning, make sure the device is turned off.
- Do not immerse the charging station (8) in water.
- Do not use strong agents, cleaning agents or cleaning products which may damage the surface.
- Wipe the charging station (8) and housing (2) with a damp cloth.
- Rinse the head (7) thoroughly. Do not compress or pull the bristles to avoid damage.

## ENVIRONMENTAL PROTECTION

- the appliance is made of materials which can be recovered or recycled
- the appliance should be handed over to the appropriate centre of collection and recycling of electrical and electronic equipment

Mode	Method of operation	Operating time
Whitening	Removes plaque, whitens front teeth and eliminates discolouration	Two minutes, impulse every 30 seconds
Cleaning	Daily plaque removal	
Gentle	Cleans sensitive teeth and gums	
Gum care	Massages gums	
Deep cleaning	Removing particularly stubborn plaque and discolouration.	Three minutes, impulse every 45 seconds

## PROBLEMS & SOLUTIONS

Problem	Possible cause	Solution
The head vibrations are too weak	The bristles are too tightly pressed against the teeth	Press the toothbrush gently against the teeth
Operating time is short	The toothbrush is not sufficiently charged: 1.it is incorrectly placed on the charging station or not in the proper spot, which causes the charging function to fail 2.the charging station is not connected to the mains.	1. place the toothbrush on the charging station properly, 2. check the connection of the charging station to the mains.
The toothbrush does not work	1.The toothbrush is not charged. 2.The toothbrush has been damaged.	1.Charge the toothbrush. 2.Contact the manufacturer's maintenance service.

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